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www.vernier.com

Vernier Graphical Analysis™ Pro

Collect and Analyze Data in Real Time for Engaged STEM Learning

Our award-winning Graphical Analysis Pro app enhances collaboration and deepens analysis of experiment data.

[Home](#) > [Products](#) > Vernier Graphical Analysis™ Pro

Inspire Your Students with Graphical Analysis Pro

With Graphical Analysis Pro, students can visualize scientific data collected from Vernier sensors—boosting their understanding and engagement using rich, real-time graphs. This easy-to-use app works on operating systems commonly found in the classroom.

Graphical Analysis Pro also includes crucial features that promote three-dimensional learning, critical thinking, student collaboration, asynchronous instruction, and deeper analysis of experiment phenomena.

Review the [comparison table](#) to see what Graphical Analysis Pro can add to your classroom.

We recognize the unique challenges you and your colleagues face during these difficult times, and we're here to help. **We are extending the free trial period for Graphical Analysis Pro from 30 days to the end of June 2022.**

[LEARN MORE](#)

[DOWNLOAD FREE VERSION](#)

[START FREE PRO TRIAL](#)

[PURCHASE NOW](#)

Data Sharing in Vernier Graphical Analysis™ Pro Overview

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www.vernier.com

Download Vernier Graphical Analysis™

With Graphical Analysis, you and your students can collect and analyze data in real time, offering an engaging learning experience.

Collect, share, and analyze sensor data with our free software for Chrome™, Windows®, macOS®, iOS, and Android™.

Ready to get started? Download the app for free.

Download for
Windows

Download for
macOS

Download on the
App Store

available in the
chrome web store

GET IT ON
Google Play


Chromebooks require the version from the Chrome Web Store (not Google Play). [Learn more »](#)

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




Solution:

Plot the points on a coordinate plane .

Calculate the means of the x -values and the y -values.

$$\bar{X} = \frac{8+2+11+6+5+4+12+9+6+1}{10} = 6.4$$

$$\bar{Y} = \frac{3+10+3+6+8+12+1+4+9+14}{10} = 7$$

Now calculate $x_i - \bar{X}$, $y_i - \bar{Y}$, $(x_i - \bar{X}) (y_i - \bar{Y})$, and $(x_i - \bar{X})^2$ for each i .

i	x_i	y_i	$x_i - \bar{X}$	$y_i - \bar{Y}$	$(x_i - \bar{X}) (y_i - \bar{Y})$	$(x_i - \bar{X})^2$
1	8	3	1.6	-4	-6.4	2.56
2	2	10	-4.4	3	-13.2	19.36
3	11	3	4.6	-4	-18.4	21.16
4	6	6	-0.4	-1	0.4	0.16
5	5	8	-1.4	1	-1.4	1.96
6	4	12	-2.4	5	-12	5.76
7	12	1	5.6	-6	-33.6	31.36
8	9	4	2.6	-3	-7.8	6.76
9	6	9	-0.4	2	-0.8	0.16
10	1	14	-5.4	7	-37.8	29.16
					$\sum_{i=1}^n (x_i - \bar{X}) (y_i - \bar{Y})$ = -131	$\sum_{i=1}^n (x_i - \bar{X})^2$ = 118.4

Calculate the slope.

$$m = \frac{\sum_{i=1}^n (x_i - \bar{X}) (y_i - \bar{Y})}{\sum_{i=1}^n (x_i - \bar{X})^2} = \frac{-131}{118.4} \approx -1.1$$

Calculate the y -intercept.

Use the formula to compute the y -intercept.


$$\begin{aligned} b &= \bar{Y} - m\bar{X} \\ &= 7 - (-1.1 \times 6.4) \\ &= 7 + 7.04 \\ &\approx 14.0 \end{aligned}$$

Use the slope and y -intercept to form the equation of the line of best fit.

The slope of the line is -1.1 and the y -intercept is 14.0 .

Therefore, the equation is $y = -1.1x + 14.0$.

Draw the line on the scatter plot.

Vaško lydymosi temperatūra



Fig 2: Graph of heating and cooling of paraffin wax

www.vernier.com

Scientific Technology for Educators & Students

For more than 40 years, we have empowered educators like you with world-class data-collection technology and innovative experiments. Your passion and dedication, along with the implementation of high-quality sensors, experiments, and resources in your classroom or laboratory, enable your students to explore science in new ways. Our mission is to provide you with the tools you need to encourage scientific curiosity in all students. See what partnering with us can do.

Elementary School

Middle School

High School

College

During these challenging times, we've got you covered.

We're supporting educators like you by extending the trial periods for our Vernier Graphical Analysis™ Pro and Vernier Video Analysis® apps to the summer, offering a free webinar for getting started with our remote learning resources, and providing free coding activities. We're here to help every step of the way!

TRY GRAPHICAL ANALYSIS PRO
FOR FREE

GET STARTED WITH
VIDEO ANALYSIS FOR FREE

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Experiment and Sample Data Library

Home > Remote Learning Solutions > Experiment and Sample Data Library

Experiment Name	Subject	Software Used for Data Analysis	Related Content
Newton's Second Law	Physics	Graphical Analysis Pro	Newton's Second Law
Newton's Third Law	Physics	Graphical Analysis Pro	Newton's Third Law
Newton's First Law	Physics	Graphical Analysis Pro	Newton's First Law
Newton's Law of Universal Gravitation	Physics	Graphical Analysis Pro	Newton's Law of Universal Gravitation
Newton's Law of Cooling	Physics	Graphical Analysis Pro	Newton's Law of Cooling
Newton's Law of Viscosity	Physics	Graphical Analysis Pro	Newton's Law of Viscosity
Newton's Law of Diffusion	Physics	Graphical Analysis Pro	Newton's Law of Diffusion
Newton's Law of Refraction	Physics	Graphical Analysis Pro	Newton's Law of Refraction
Newton's Law of Reflection	Physics	Graphical Analysis Pro	Newton's Law of Reflection
Newton's Law of Diffraction	Physics	Graphical Analysis Pro	Newton's Law of Diffraction
Newton's Law of Interference	Physics	Graphical Analysis Pro	Newton's Law of Interference
Newton's Law of Polarization	Physics	Graphical Analysis Pro	Newton's Law of Polarization
Newton's Law of Dispersion	Physics	Graphical Analysis Pro	Newton's Law of Dispersion
Newton's Law of Absorption	Physics	Graphical Analysis Pro	Newton's Law of Absorption
Newton's Law of Emission	Physics	Graphical Analysis Pro	Newton's Law of Emission
Newton's Law of Scattering	Physics	Graphical Analysis Pro	Newton's Law of Scattering
Newton's Law of Diffusion	Physics	Graphical Analysis Pro	Newton's Law of Diffusion
Newton's Law of Refraction	Physics	Graphical Analysis Pro	Newton's Law of Refraction
Newton's Law of Reflection	Physics	Graphical Analysis Pro	Newton's Law of Reflection
Newton's Law of Diffraction	Physics	Graphical Analysis Pro	Newton's Law of Diffraction
Newton's Law of Interference	Physics	Graphical Analysis Pro	Newton's Law of Interference
Newton's Law of Polarization	Physics	Graphical Analysis Pro	Newton's Law of Polarization
Newton's Law of Dispersion	Physics	Graphical Analysis Pro	Newton's Law of Dispersion
Newton's Law of Absorption	Physics	Graphical Analysis Pro	Newton's Law of Absorption
Newton's Law of Emission	Physics	Graphical Analysis Pro	Newton's Law of Emission
Newton's Law of Scattering	Physics	Graphical Analysis Pro	Newton's Law of Scattering

Bring Science to Your Students at Home

To assist you as you move from the classroom or laboratory setting to online teaching, we are offering the Vernier Experiment and Sample Data Library—a free collection of over 300 experiments with sample data covering many subjects that you can distribute to your students at home.


Students will use the experiment data we've collected to perform their own analysis.

The Vernier Experiment and Sample Data Library is available through June 30, 2022. We have a growing library of experiments with data *and* videos included with Graphical Analysis Pro (free trial available) to bring science to your students at home.

[BROWSE LIBRARY](#)

How do I use the Experiment and Sample Data Library?

Watch our overview tutorial to learn how to use our experiments and data with your students.

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www.vernier.com

Experiment and Sample Data Library Index

Biological Sciences | Chemistry | Physics | Elementary/Middle/Physical Science

Rev. 7.23.2020

[Download Instructions for Using the Library Here](#)

Index of Biological Sciences	Software Used for Data Analysis	Related videos (where available)
Follow the links below to access the experiments. You may download all the biology, human physiology, and agricultural science experiments, a specific lab book, or individual experiments. Each individual experiment folder will contain the experiment instructions as a Word doc and one or more sample data files.	Use the software indicated here to view and analyze the experiment's sample data. Some titles are free, others require licenses. If your school or department doesn't own a license, a free demo is available.	Note that while the science and the experimental setup shown in these videos match the written instructions, the Vernier equipment and software may be older and may not match. Students should follow the written instructions provided whenever possible.
ALL BIOLOGICAL SCIENCES		
Advanced Biology with Vernier (BIO-A)		
BIO-A 01A Membrane Diffusion	Logger Pro 3	
BIO-A 02 Enzyme (O2)	Logger Pro 3	https://youtu.be/2fNu0ARUpDQ
BIO-A 02 Enzyme (Pressure)	Logger Pro 3	https://youtu.be/UtFzkQyrAIU
BIO-A 04B Photosynthesis	Logger Pro 3	
BIO-A 05 Cell Respiration (CO2)	Logger Pro 3	https://youtu.be/xQ-kEV7UCno
BIO-A 06B DNA Fingerprinting	Logger Pro 3	
BIO-A 09 Transpiration	Logger Pro 3	https://youtu.be/oeah4AMH5sU
BIO-A 10A Blood Pressure	Logger Pro 3	https://youtu.be/wbl64nAUmlw
BIO-A 10B Heart Rate & Fitness	Logger Pro 3	https://youtu.be/J7_R0_wUPkk
BIO-A 12A Dissolved Oxygen	Logger Pro 3	https://youtu.be/zHHRmqWuTDo
BIO-A 12B Primary Productivity	Logger Pro 3	
BIO-A 13 Visible Spectra of Plant Pigments	Logger Pro 3	https://youtu.be/UxZfrcdOCK8
BIO-A 14 Determination of Chlorophyll in Olive Oil	Logger Pro 3	
BIO-A 15 Enzyme Analysis using Tyrosinase	Logger Pro 3	
BIO-A 16 Intro to Neurotransmitters using AChE	Logger Pro 3	
BIO-A 17 Macromolecules Experiments with Protein	Logger Pro 3	
Biology with Vernier (BWV)		
BWV 01 Energy in Food	Graphical Analysis 4	

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Graphical Analysis 1

Energy in Food

Food supplies energy for all animals – without it we could not live. The quantity of energy stored in food is of great interest to humans. The energy your body needs for running, sitting, and thinking comes from the food you eat. Not all foods contain the same amount of energy, nor are all foods equally nutritious. An average person should consume a minimum of 2,000 kilocalories per day, which is equivalent to 8,400 kilojoules. Calories and joules are both units of energy. We will use joules in this lab since it is the accepted SI metric standard.

You can determine energy content of food by burning a portion of it and capturing the heat released in a known amount of water. This technique is called calorimetry. The energy content of the food is the amount of heat produced by the combustion of 1 gram of the food substance (i.e., measured in kilojoules per gram (kJ/g)).

OBJECTIVES

- Use a Temperature Probe to measure temperature changes.
- Monitor the energy gained off by food as it burns.
- Determine and compare the energy content of different foods.

Figure 1

Energy in Food

MATERIALS

Chrombook, computer, or mobile device
Graphical Analysis 4 app
Go Direct Temperature
100 mL graduated cylinder
balance
food holder
two food samples
matches
ring stand
10 cm ring support
small can
two stirring rods
utility clamp
cold water
wooden splint
graphing

PROCEDURE

- Obtain and wear goggles.
- Launch Graphical Analysis. Connect the Temperature Probe to your Chrombook, computer, or mobile device.
- Set up the Data collection mode.
 - Click or tap Mode to open Data Collection Settings.
 - Change Rate to 0.5 samples.
 - Set Data Collection to 600 s. Click or tap Done.
- Obtain a sample of food and a food holder. Mount the food onto the food holder so that it can burn without damaging the holder (see Figure 1). Determine and record the initial mass of the food sample and food holder in Table 1. **Caution:** Do not eat or drink in the laboratory.
- Setup the apparatus shown in Figure 1.
 - Determine the mass of an empty can. Record the value in Table 1.
 - Place 90 mL of cold water into the can.
 - Determine and record the mass of the can plus the water.
 - Insert a stirring rod through the holes in the top of the can. Use a 10 cm ring and ring stand to suspend the can about 2.5 cm (1 inch) above the food sample.
 - Use a utility clamp to suspend the Temperature Probe in the water as shown in Figure 1. The probe should not touch the bottom or side of the can.
- Click or tap Collect to start data collection.

Energy in Food

DATA

Table 1		
	Sample 1	Sample 2
Point used:		
Mass of empty can (g)		
Mass of can plus water (g)		
Minimum temperature of water (°C)		
Maximum temperature of water (°C)		
Initial mass of food and holder (g)		
Final mass of food and holder (g)		

Table 2		
	Sample 1	Sample 2
Mass of water (g)		
ΔT of water (°C)		
Change in mass of food (g)		
Energy gained by water (J)		
Energy content of food (kJ/g)		

Energy in Food

Table 3: Class Results

Lab group	Sample 1	Sample 2
	Energy content (kJ/g)	Energy content (kJ/g)
Average		

PROCESSING THE DATA

Record the following calculations in Table 2.

- Calculate the mass of water for each sample. Show your calculations.
- Calculate the change in mass of each food sample burned. Show your calculations.
- Calculate the change in the temperature of the water, ΔT. Show your calculations.
- Calculate the energy gained by the heated water using the equation below. Show your calculations.

$$\text{Energy gained by water} = (\text{mass of water}) \times (\Delta T \text{ of water}) \times (4.18 \text{ J/g}^\circ\text{C})$$
- Convert the energy you calculated in Step 4 to kilojoules (1 kJ = 1000 J). Show your calculation.
- Use your answer in Step 2 and 5 to calculate the energy content of each food sample (in kJ/g).

$$\text{Energy content of food} = \text{Energy gained by water} / \text{change of mass of food}$$
- Record your results and the results of other groups in Table 3.

Energy in Food

- Use a match to light the food sample. Quickly position the burning food sample directly below the center of the water-filled can. **Caution:** Always keep hair and clothing away from open flames.
- Stir the water slowly and continuously using the stirring rod while monitoring the temperature. If the temperature of the water exceeds 60°C, blow out the flame.
- If the food is still burning after 4 minutes, blow out the flame. Do not stop data collection yet.
- Once the water temperature begins to decrease, stop data collection.
- After data collection has stopped, analyze the graph to determine the maximum and minimum temperatures of the water.
 - Click or tap Graph Tools, **MC**, and choose View Statistics.
 - Record the Maximum (Final) and Minimum (Starting) temperature values recorded during data collection in your data table.
 - Dismiss the Statistics box.
- Determine the final mass of the food and holder. Record the mass in Table 1.
- Clean off the food holder and empty the can of water. Discard all burnt matches and food samples as directed by your instructor.
- Repeat Steps 4–11 for the second food sample. **Note:** The previous data set is automatically saved.
- A good way to compare the two samples is to view both sets of data on one graph.
 - To do this, click or tap the y-axis label and select the data sets you want to display. Dismiss the box to view the graph.
 - Examine the data points along the displayed curve. Click or tap the graph to maximize the data. **Note:** You can also adjust the Examine line by dragging the line.


QUESTIONS

- Which of the foods has the greatest energy content?
- Which of the tested foods is the best energy source? Why?
- What was the original energy source of the foods tested?
- Why might some foods with a lower energy content be better energy sources than other foods with a higher energy content?
- Would you expect the energy content values that you measured to be close to the value listed in dietary books? Why?

EXTENSION

Determine the energy content of other combustible foods.

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